

Product Datasheet

Ubiquilin1 Antibody Blocking peptide (orb1450994)

Description	Ubiquilin1 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1939522 was selected from the Center region of human Ubiquilin1 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q9UMX0
Expiration Date	12 months from date of receipt.