

Product Datasheet

UBE2V1 Antibody Blocking peptide (orb1451030)

| | |
|--------------------------|---|
| Description | UBE2V1 Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence was selected from the N-term region of human UBE2V1 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | Q13404 |
| Expiration Date | 12 months from date of receipt. |