

Product Datasheet

USP2 (Ctr S260) Antibody Blocking peptide (orb1451084)

Description	USP2 (Ctr S260) Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931704 was selected from the Ctr region of human USP2 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	O75604
Expiration Date	12 months from date of receipt.