

Product Datasheet

USP2 (L523) Antibody Blocking peptide (orb1451085)

Description USP2 (L523) Antibody Blocking peptide

Form/Appearance Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Note For research use only

Application notes Synthetic peptide

Protein Sequence The synthetic peptide sequence was selected from the C-term region of human USP2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Uniprot ID [O75604](#)

Expiration Date 12 months from date of receipt.