

## Product Datasheet

### UBE2L6 Antibody Blocking peptide (orb1451115)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | UBE2L6 Antibody Blocking peptide   |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.   |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Synthetic peptide  |
| <b>Protein Sequence</b>  | The synthetic peptide sequence used to generate the antibody orb1931722 was selected from the N-term region of human UBE2L6 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>Uniprot ID</b>        | <b>O14933</b>  |
| <b>Expiration Date</b>   | 12 months from date of receipt.  |