

Product Datasheet

HIP2 Antibody Blocking peptide (orb1451122)

Description	HIP2 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1931726 was selected from the Center region of human HIP2 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	P61086
Expiration Date	12 months from date of receipt.