

## Product Datasheet

# USP24 Antibody (Y1427) Blocking peptide (orb1451923)

<b>Description</b>	USP24 Antibody (Y1427) Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence was selected from the C-term region of human USP24. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
<b>Uniprot ID</b>	<b>Q9UPU5</b>
<b>Expiration Date</b>	12 months from date of receipt.