



Product Datasheet

Zap-70 (Ab-493) Antibody (orb14568)

Description	Zap-70 (Ab-493) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P43403
Expiration Date	12 months from date of receipt.



Western blot analysis of Zap-70(Ab-493) ...



Immunohistochemical analysis of formalin...