

## Product Datasheet

### MCOLN1 Antibody Blocking peptide (orb1456990)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | MCOLN1 Antibody Blocking peptide   |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.   |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | Synthetic peptide  |
| <b>Protein Sequence</b>  | The synthetic peptide sequence was selected from the C-term region of MCOLN1. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>Uniprot ID</b>        | <a href="#">Q9GZU1</a>   |
| <b>Expiration Date</b>   | 12 months from date of receipt.  |