



Product Datasheet

STAT1 (Ab-727) Antibody (orb14572)

Description	STAT1 (Ab-727) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P42224
Expiration Date	12 months from date of receipt.



Western blot analysis of STAT1(Ab-727) a...



Immunohistochemical analysis of formalin...