

Product Datasheet

Bmf BH3 Domain Antibody Blocking peptide (orb1457487)

Description	Bmf BH3 Domain Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1936936 was selected from the region of human Bmf BH3 Domain. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q96LC9
Expiration Date	12 months from date of receipt.