

Product Datasheet

BCL2A1 (A1) (BH3 Domain Specific) Antibody Blocking peptide (orb1457588)

| | |
|--------------------------|---|
| Description | BCL2A1 (A1) (BH3 Domain Specific) Antibody Blocking peptide |
| Form/Appearance | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| Note | For research use only |
| Application notes | Synthetic peptide |
| Protein Sequence | The synthetic peptide sequence used to generate the antibody orb1937011 was selected from the region of human A1 BH3 Domain. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| Uniprot ID | Q16548 |
| Expiration Date | 12 months from date of receipt. |