



Product Datasheet

AMPK α 1/AMPK α 2 (Ab183/172) Antibody (orb14583)


Description	AMPK α 1/AMPK α 2 (Ab183/172) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q13131
Expiration Date	12 months from date of receipt.



Western blot analysis of AMPK α 1/AMPK α 2(A...



Confocal immunofluorescent analysis of H...



Immunohistochemical analysis of formalin...