

## Product Datasheet

### IL1F5 Antibody Blocking peptide (orb1458336)

|                        |  |
|------------------------|--|
| <b>Description</b>     | IL1F5 Antibody Blocking peptide  |
| <b>Form/Appearance</b> | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed. |
| <b>Note</b>            | For research use only  |
| <b>Uniprot ID</b>      | <a href="#">Q9UBH0</a>   |
| <b>Expiration Date</b> | 12 months from date of receipt.  |