

## Product Datasheet

### SET2 Antibody Blocking peptide (orb1458709)

<b>Description</b>	SET2 Antibody Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence used to generate the antibody orb1937827 was selected from the N-term region of human SET2. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
<b>Uniprot ID</b>	<a href="#">Q9BYW2</a>
<b>Expiration Date</b>	12 months from date of receipt.