



Product Datasheet

APP (Ab-668) Antibody (orb14595)

Description	APP (Ab-668) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05067
Expiration Date	12 months from date of receipt.



Western blot analysis of APP (Ab-668) ant...



Confocal immunofluorescent analysis of M...