

## Product Datasheet

# SAE2 (UBA2) Antibody (A600) Blocking peptide (orb1460110)

<b>Description</b>	SAE2 (UBA2) Antibody (A600) Blocking peptide
<b>Form/Appearance</b>	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
<b>Note</b>	For research use only
<b>Application notes</b>	Synthetic peptide
<b>Protein Sequence</b>	The synthetic peptide sequence was selected from the C-term region of human UBA2 . A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
<b>Uniprot ID</b>	<b>Q9UBT2</b>
<b>Expiration Date</b>	12 months from date of receipt.

---

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)