

## Product Datasheet

### UBE1C Antibody Blocking peptide (orb1460136)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | UBE1C Antibody Blocking peptide   |
| <b>Form/Appearance</b>   | Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.  |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Synthetic peptide   |
| <b>Protein Sequence</b>  | The synthetic peptide sequence was selected from the N-term region of human UBE1C. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay. |
| <b>Uniprot ID</b>        | <a href="#">Q8TBC4</a>  |
| <b>Expiration Date</b>   | 12 months from date of receipt.   |