

Product Datasheet

Me2-H3(K9) Antibody Blocking peptide (orb1460261)

Description Me2-H3(K9) Antibody Blocking peptide

Form/Appearance Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.

Note For research use only

Application notes Synthetic peptide

Protein Sequence The synthetic peptide sequence was selected from the 2Me region of human Methyl-H3-mK9(2Me). A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.

Uniprot ID [Q6NXT2](#)

Expiration Date 12 months from date of receipt.