

Product Datasheet

JHDM1a Antibody Blocking peptide (orb1460347)

Description	JHDM1a Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence used to generate the antibody orb1938986 was selected from the Center region of human JHDM1a/FBXL11. A 10 to 100 fold molar excess to antibody is recommended. Precise conditions should be optimized for a particular assay.
Uniprot ID	Q9Y2K7
Expiration Date	12 months from date of receipt.