

Product Datasheet

KY02111 (orb146070)

| | |
|------------------------|---|
| Description | KY02111 promotes differentiation of hPSCs to cardiomyocytes by inhibiting WNT signaling. |
| Note | For research use only |
| Formula | $C_{18}H_{17}CLN_2O_3S$ |
| MW | 376.8572 |
| CAS Number | [1118807-13-8] |
| SMILES | <chem>C1C([H])=C([H])C2=C(C=1[H])SC(=N2)N([H])C(C([H])([H])C([H])([H])C1C([H])=C([H])C=C(C=1[H])OC([H])([H])[H])OC([H])([H])[H])=O</chem> |
| Expiration Date | 12 months from date of receipt. |