

Product Datasheet

ZSTK474 (orb146082)

Description	ZSTK474 is a potent pan-class I PI3K inhibitor with IC50 of 37 nM.
Note	For research use only
Formula	$C_{19}H_{21}F_2N_7O_2$
MW	417.41
CAS Number	[475110-96-4]
SMILES	<chem>FC(C1=NC2=C(N1C3=NC(N4CCOCC4)=NC(N5CCOCC5)=N3)C=CC=C2)F</chem>
Expiration Date	12 months from date of receipt.