

Product Datasheet

PRMT1 Antibody Blocking peptide (orb1460834)

Description	PRMT1 Antibody Blocking peptide
Form/Appearance	Peptides are lyophilized in a solid powder format. Peptides can be reconstituted in solution using the appropriate buffer as needed.
Note	For research use only
Application notes	Synthetic peptide
Protein Sequence	The synthetic peptide sequence is selected from aa 23-40 of human PRMT1.
Uniprot ID	Q99873
Expiration Date	12 months from date of receipt.