

## Product Datasheet

### PCI-34051 (orb146092)

<b>Description</b>	PCI-34051 is a potent and specific HDAC8 inhibitor with IC50 of 10 nM.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{17}H_{16}N_2O_3$
<b>MW</b>	296.32
<b>CAS Number</b>	[950762-95-5]
<b>SMILES</b>	<chem>COC1=CC=C(C=C1)CN2C=CC3=C2C=C(C=C3)C(=O)NO</chem>
<b>Expiration Date</b>	12 months from date of receipt.