

Product Datasheet

Caldesmon, HMW (h-Caldesmon) antibody (orb1460961)

Description	Rabbit monoclonal antibody to Caldesmon, HMW (h-Caldesmon)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q05682
Expiration Date	12 months from date of receipt.