

Product Datasheet

HLA-DP antibody (orb1461244)

Description	Mouse monoclonal antibody to HLA-DP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P04440
Expiration Date	12 months from date of receipt.