

## Product Datasheet

### Dovitinib (TKI258) Lactate (orb146142)

<b>Description</b>	Dovitinib (TKI258, CHIR258) is a novel multi-target inhibitor for Flt3, c-Kit, FGFR1/3, VEGFR1/2/3, ...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{21}H_{21}N_6O_4 \cdot C_3H_6O_3 \cdot H_2O$
<b>Purity</b>	> 98%, Standard References
<b>MW</b>	500.52266
<b>CAS Number</b>	[915769-50-5]
<b>SMILES</b>	<chem>O=C1NC2=C(C(F)=CC=C2)C(N)=C1C3=NC4=CC=C(N5CCN(C)CC5)C=C4N3.CC(O)C(O)=O.[H]O[H]</chem>
<b>Expiration Date</b>	12 months from date of receipt.