

## Product Datasheet

### Darifenacin (orb146146)

<b>Description</b>	Darifenacin is a selective M3 muscarinic receptor antagonist with pKi of 8.9.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{28}H_{30}N_2O_2$
<b>MW</b>	426.5
<b>CAS Number</b>	[133099-04-4]
<b>SMILES</b>	<chem>NC(C(C1=CC=CC=C1)(C2=CC=CC=C2)[C@H](CC3)CN3CCC4=CC(CCO5)=C5C=C4)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.