



Product Datasheet

S6 Ribosomal Protein (Ab-235) Antibody (orb14615)


Description	S6 Ribosomal Protein (Ab-235) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P62753
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of M...



Western blot analysis of S6 Ribosomal Pr...



Immunohistochemical analysis of paraffin...