

Product Datasheet

Nolatrexed (AG-337) (orb146163)

Description	Nolatrexed (AG-337), inhibits purified recombinant human thymidylate synthase (TS) with a K_i of 11 n...
Note	For research use only
Formula	$C_{14}H_{14}Cl_2N_4OS$
MW	357.25
CAS Number	[152946-68-4]
SMILES	<chem>CC1=C(SC2=CC=NC=C2)C(C3=O)=C(N=C(N)N3)C=C1</chem>
Expiration Date	12 months from date of receipt.