

## Product Datasheet

# GRP94 / HSP90B1 (Endoplasmic Reticulum Marker) - With BSA and Azide antibody (orb1461866)

<b>Description</b>	Rat monoclonal antibody to GRP94 / HSP90B1 (Endoplasmic Reticulum Marker) - With BSA and Azide
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Bovine, Canine, Equine, Guinea pig, Hamster, Human, Mouse, Other, Porcine, Rabbit, Rat, Sheep
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC-P, WB
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P14625</b>
<b>Expiration Date</b>	12 months from date of receipt.