

Product Datasheet

CD71 / Transferrin Receptor (TFRC) - With BSA and Azide antibody (orb1461900)

Description	Mouse monoclonal antibody to CD71 / Transferrin Receptor (TFRC) - With BSA and Azide
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P02786
Expiration Date	12 months from date of receipt.