

Product Datasheet

LMK-235 (orb146208)

Description	LMK-235 is a selective histone deacetylase (HDAC) 4 and HDAC5 inhibitor (IC50 values are 4.22, 11.9,...)
Note	For research use only
Formula	$C_{15}H_{22}N_2O_4$
MW	294.3462
CAS Number	[1418033-25-6]
SMILES	<chem>O(C([H])([H])C([H])([H])C([H])([H])C([H])([H])C([H])([H])C([H])([H])C(N([H])O[H])=O)N([H])C1C([H])=C(C([H])([H])[H])C([H])=C(C([H])([H])[H])C=1[H])=O</chem>
Expiration Date	12 months from date of receipt.