

Product Datasheet

Neurofilament, phospho (NF-H) (Neuronal Marker) - With BSA and Azide antibody (orb1462121)

Description	Mouse monoclonal antibody to Neurofilament, phospho (NF-H) (Neuronal Marker) - With BSA and Azide
Species/Host	Mouse
Reactivity	Bovine, Feline, Guinea pig, Human, Mouse, Porcine, Rabbit, Rat, Rodent
Conjugation	Unconjugated
Tested Applications	FC, IHC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P12036
Expiration Date	12 months from date of receipt.