

## Product Datasheet

### Enzastaurin (LY317615) (orb146243)

<b>Description</b>	Enzastaurin (LY317615) is a potent PKC $\beta$ selective inhibitor with IC50 of 6 nM, 6- to 20-fold select...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>32</sub> H <sub>29</sub> N <sub>5</sub> O <sub>2</sub>
<b>MW</b>	515.60496
<b>CAS Number</b>	[170364-57-5]
<b>SMILES</b>	<chem>O=C(C(C1=CN(C)C2=C1C=CC=C2)=C3C4=CN(C5CCN(CC6=NC=CC=C6)CC5)C7=C4C=CC=C7)NC3=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.