

## Product Datasheet

# HLA-DRB (MHC II) - With BSA and Azide antibody (orb1462444)

<b>Description</b>	Mouse monoclonal antibody to HLA-DRB (MHC II) - With BSA and Azide
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Canine, Human, Primate, Rodent
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P01911</a>
<b>Expiration Date</b>	12 months from date of receipt.