

Product Datasheet

(+)PD 128907 (orb146245)

Description	(+)PD 128907 is a potent D3 dopamine receptor agonist (K _i = 2.3 nM). Displays 18 - 200-fold selectiv...
Note	For research use only
Formula	C ₁₄ H ₁₉ NO ₃ ·HCl
Purity	> 95%
MW	285.77
CAS Number	[300576-59-4]
SMILES	OC1=CC=C2C([C@@]3([H])OCCN(CCC)[C@]3([H])CO2)=C1.[H]Cl
Expiration Date	12 months from date of receipt.