

Product Datasheet

HLA-DRB (MHC II) - With BSA and Azide antibody (orb1462450)

Description	Mouse monoclonal antibody to HLA-DRB (MHC II) - With BSA and Azide
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P01911
Expiration Date	12 months from date of receipt.