

Product Datasheet

Ceritinib(LDK378) (orb146278)

Description	LDK378 is a highly selective ALK inhibitor with IC50 of 2 nM
Note	For research use only
Formula	$C_{28}H_{36}N_5O_3SCL$
MW	558.13514
CAS Number	[1032900-25-6]
SMILES	<chem>CC(C)OC1=CC(C2CCNCC2)=C(C)C=C1NC3=NC=C(Cl)C(NC4=CC=CC=C4S(=O)(C(C)C)=O)=N3</chem>
Expiration Date	12 months from date of receipt.