



Product Datasheet

NFκB-p105 (Ab-932) Antibody (orb14629)

Description	NFκB-p105 (Ab-932) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P19838
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of NFκB-p105(Ab-932) in HeLa cells. The blot shows a single band at approximately 100 kD. The lanes are labeled 'Hela' and 'Peptide - +'. The band intensity is significantly reduced in the '+' lane compared to the '-' lane.