

## Product Datasheet

# TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) - With BSA and Azide antibody (orb1463029)

<b>Description</b>	Mouse monoclonal antibody to TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) - With BSA and Azide
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Bovine, Canine, Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC, IF, IHC-P
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Expiration Date</b>	12 months from date of receipt.