

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) - With BSA and Azide antibody (orb1463031)

Description	Mouse monoclonal antibody to TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) - With BSA and Azide
Species/Host	Mouse
Reactivity	Bovine, Canine, Hamster, Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.