

Product Datasheet

Arginase-1 antibody (orb1463153)

Description	Rabbit monoclonal antibody to Arginase-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Entrez	383
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.