

Product Datasheet

Galectin-3 antibody (orb1463207)

Description	Rabbit monoclonal antibody to Galectin-3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Target	LGALS3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	26152
Uniprot ID	P17931
Entrez	3958
Dilution Range	IHC: Ready-to-use
Expiration Date	12 months from date of receipt.