

---

## Product Datasheet

### Membrane Protein (SARS-CoV-2) antibody (orb1463313)

**Description**

Membrane protein (E) is one of the

glycoprotein that is the m...

**Species/Host**

Goat

**Reactivity**

Virus

**Conjugation**

Unconjugated

**Tested Applications**

WB

**Immunogen**

Affinity purified recombinant fusion protein using the C-terminal of Membrane protein (residues 100 to stop) and produced in E. coli.

**Target**

Membrane Protein (SARS-CoV-2)

**Preservatives**

PBS, 20% glycerol and 0.05% sodium azide

**Form/Appearance**

Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

**Concentration**

3 mg/ml

**Storage**

Store at -20°C for long-term storage.  
Store at 2-8°C for up to one month.

**Note**

For research use only

**Application notes**

The antibody solution should be gently mixed before use.

**Isotype**

IgG

**Clonality**

Polyclonal

**Dilution Range**

WB:1:500-1:2,000

**Expiration Date**

12 months from date of receipt.