

Product Datasheet

NSP3 (SARS-CoV-2) antibody (orb1463325)

Description

Nsp3 is part of the multifunctional protein

in the transc...

Species/Host

Goat

Reactivity

Virus

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Affinity purified recombinant fusion protein using the C-terminal of Nsp3 (residues 1846 to stop) and produced in E. coli.

Target

NSP3 (SARS-CoV-2)

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

3 mg/ml

Storage

Store at -20°C for long-term storage.
Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB:1:500-1:2,000

Expiration Date

12 months from date of receipt.