

Product Datasheet

NSP1 (SARS-CoV-2) antibody (orb1463327)

Description	Host translation inhibitor Nsp1 is part of the multifunctional protein replicase polyprotein 1ab tha...
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Affinity purified recombinant fusion protein Nsp1 of SARS-CoV-2 produced in E. coli.
Target	NSP1 (SARS-CoV-2)
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:2,000
Expiration Date	12 months from date of receipt.