
Product Datasheet

Spike RBD Domain (SARS-CoV-2) antibody (orb1463332)

Description

The spike or S glycoprotein is a

molecular weight of about 150 kDa
foun...

Species/Host

Goat

Reactivity

Virus

Conjugation

Unconjugated

**Tested
Applications**

WB

Immunogen

Affinity purified recombinant fusion protein using the spike protein RBD domain (residues 320 to 420) and produced in E. coli.

Target

Spike RBD Domain (SARS-CoV-2)

Preservatives

PBS, 20% glycerol and 0.05% sodium azide

Form/Appearance

Polyclonal antibody supplied as a 100 µl (3 mg/ml) aliquot in PBS, 20% glycerol and 0.05% sodium azide. This antibody is epitope-affinity purified from goat antiserum.

Concentration

3 mg/ml

Storage

Store at -20°C for long-term storage.
Store at 2-8°C for up to one month.

Note

For research use only

Application notes

The antibody solution should be gently mixed before use.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

WB:1:500-1:2,000

Expiration Date

12 months from date of receipt.