

Product Datasheet

Streptavidin antibody (orb1463338)

Description	Goat polyclonal antibody to streptavidin. Streptavidin is a tetrameric bacterial protein isolated fr...
Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide derived from within residues 50 aa to the C-terminus of streptavidin protein.
Target	Streptavidin
Preservatives	PBS, 20% glycerol and 0.05% sodium azide
Form/Appearance	Polyclonal antibody supplied as a 200 µl (2 mg/ml) aliquot in NaHCO ₃ 100 mM pH 9. This antibody is epitope-affinity purified from goat antiserum and does not contain preservatives.
Concentration	3 mg/ml
Storage	Store at -20°C for long-term storage. Store at 2-8°C for up to one month.
Note	For research use only
Application notes	The antibody solution should be gently mixed before use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB:1:500-1:5,000
Expiration Date	12 months from date of receipt.