

Product Datasheet

HA antibody (orb1463346)

Description	Human monoclonal antibody to HA
Species/Host	Human
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	BLI, ELISA
Immunogen	This antibody was raised by cloning antibody sequences from B cells isolated from people immunized with the H5 DNA vaccine.
Target	HA
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The affinity of this antibody to H1/1999 HA was assessed. This was done by using the human version of this antibody for Bio-Layer Interferometry on H1/1999 HA. Further, an ELISA was performed using the human version of this antibody on HA from influenza (Whittle et al, 2014; pmid:24501410).
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	A5A5M5
Expiration Date	12 months from date of receipt.